

565 Nucla Way, Unit B
Aurora, Colorado 80011-9319
1-800-525-2130
Phone 303-366-1804
Fax 303-367-5118
www.doweaver.com

Product Usage After Opening

All Weaver and Company products have a three (3) year shelf-life from the date of manufacture. This shelf-life is not diminished in any way after the product container has been opened, as long as proper storage of the container takes place.


We recommend the following storage conditions of our products:

- Keep container tightly closed when not in use;
- Store container at room temperature.

Nuprep® Skin Prep Gel and Ten20® Conductive Paste are not intended for single use. They are also not sterile devices. This indicates that once a tube or jar has been opened and used on a patient for the first time, the tube or jar can then be closed and kept until ready to be used on the next patient. The three (3) year shelf-life indicates that the product is usable during that entire three-year span, even while the container may be opened and closed numerous times; the product does not need to be used within a certain period of time after being opened for the first time.

If you have any further questions or concerns about the usability of our products, please do not hesitate to contact us.

Thank you for your continued support of Weaver and Company products.



Chris Cooper
Chief Executive Officer
Weaver and Company

11/2/2012
Date